

TECHNICAL DATA SHEET

<u>PRODUCT</u>	Crosstex Duo-Check Sterilization Pouches																
<u>MATERIALS</u>	<p>Paper – 60gsm Medical Grade Paper (55%)</p> <p>Plastic – 2mL Polyester Cast Polypropylene (45%) (Blue tint)</p> <p>Adhesive Tape – ½ mil polyester film coated sides with 1.4 mils of high tack, shear acrylic adhesive that exhibits good temperature resistance.</p> <p>Indicators - Internal/External printed directly on pouch</p>																
<u>PRODUCT/ PACKAGING REQUIREMENTS</u>	<p>Duo-Check Sterilization Self-Sealing Pouches have an extra wide chevron seal that protects against instrument breaching. The pouch is made of medical grade paper with a transparent plastic film that facilitates both easy instrument viewing and identification of punctures or tears in film.</p> <p>The wide self-seal adhesive strip easily and securely seals pouch, providing reliable seal until pouch is opened. The Perforated Fold facilitates easy, accurate folding of adhesive strip to form airtight seal until pouch is opened.</p> <p>Duo-Check Sterilization Self-Sealing Pouches are single use pouches. High grade pouch material is designed to allow sterilization agents (steam, EtO) to penetrate inside during the sterilization cycle. After one use, pores in the pouch close to maintain sterility until pouch is opened for use.</p> <p>It is suggested to insert into the pouch an internal indicator strip to confirm sterility inside. Indicators change color with sterilization heat and are easily visible however, they reflect only that the instruments have been processed and not sterilized.</p> <table><tr><td>Item # SCZ - 2" x 8" – 200/Box</td><td>Item # SCB1X – 4.25" x 11" - 200/Box</td></tr><tr><td>Item # SCXX – 2.25" x 4" – 200/Box</td><td>Item # SCM – 5.25" x 10" - 200/Box</td></tr><tr><td>Item # SCX – 2.75" x 9" - 200/Box</td><td>Item # SCL515 – 5.25" x 15" - 200/Box</td></tr><tr><td>Item # SCQ – 3.25" x 9" - 200/Box</td><td>Item # SCL – 7.5" x 13" - 200/Box</td></tr><tr><td>Item # SCXS – 3.5" x 5.25" - 200/Box</td><td>Item # SCL816 – 8" x 16" - 200/Box</td></tr><tr><td>Item # SCS – 3.5" x 9" - 200/Box</td><td>Item # SCL1015 - 10" x 15" - 200/Box</td></tr><tr><td>Item # SCS5 – 3.5" x 9" - 500/Box</td><td>Item # SCL1215 - 12" x 15" - 200/Box</td></tr><tr><td>Item # SCL322 – 3.5" x 22" - 100/Box</td><td>Item # SCL1218 - 12" x 18" - 200/Box</td></tr></table>	Item # SCZ - 2" x 8" – 200/Box	Item # SCB1X – 4.25" x 11" - 200/Box	Item # SCXX – 2.25" x 4" – 200/Box	Item # SCM – 5.25" x 10" - 200/Box	Item # SCX – 2.75" x 9" - 200/Box	Item # SCL515 – 5.25" x 15" - 200/Box	Item # SCQ – 3.25" x 9" - 200/Box	Item # SCL – 7.5" x 13" - 200/Box	Item # SCXS – 3.5" x 5.25" - 200/Box	Item # SCL816 – 8" x 16" - 200/Box	Item # SCS – 3.5" x 9" - 200/Box	Item # SCL1015 - 10" x 15" - 200/Box	Item # SCS5 – 3.5" x 9" - 500/Box	Item # SCL1215 - 12" x 15" - 200/Box	Item # SCL322 – 3.5" x 22" - 100/Box	Item # SCL1218 - 12" x 18" - 200/Box
Item # SCZ - 2" x 8" – 200/Box	Item # SCB1X – 4.25" x 11" - 200/Box																
Item # SCXX – 2.25" x 4" – 200/Box	Item # SCM – 5.25" x 10" - 200/Box																
Item # SCX – 2.75" x 9" - 200/Box	Item # SCL515 – 5.25" x 15" - 200/Box																
Item # SCQ – 3.25" x 9" - 200/Box	Item # SCL – 7.5" x 13" - 200/Box																
Item # SCXS – 3.5" x 5.25" - 200/Box	Item # SCL816 – 8" x 16" - 200/Box																
Item # SCS – 3.5" x 9" - 200/Box	Item # SCL1015 - 10" x 15" - 200/Box																
Item # SCS5 – 3.5" x 9" - 500/Box	Item # SCL1215 - 12" x 15" - 200/Box																
Item # SCL322 – 3.5" x 22" - 100/Box	Item # SCL1218 - 12" x 18" - 200/Box																
<u>STORAGE REQUIREMENTS</u>	Store in a cool, dry place																
<u>SHELF LIFE</u>	When properly sealed, sterility is maintained up to (3) three years provided the integrity of the pouch has not been compromised; i. e. pouch becomes wet, damaged or open during storage. Caution – Wet pouches may not preserve sterility of contents.																
<u>SINGLE-USE POUCHES</u>	Sterilization pouches are single use only items. Do not re-use.																

Duo-Check Sterilization Self-Sealing Pouches are not considered hazardous waste and can be disposed of in regular trash (40 CFR 261.24 tested using the TCLP, EPA Publication SW-846)