

## TECHNICAL DATA SHEET

<b><u>PRODUCT</u></b>	Crosstex CSR Sterilization Wrap
<b><u>MATERIALS</u></b>	<ul style="list-style-type: none"> <li>• 52 gsm cellulose fibers reinforced with synthetic binders</li> <li>• Cellulose pulp &gt;80%*</li> <li>• Binders &lt;20%</li> <li>• Sizing Agent 0.3%</li> <li>• Wet Strength Agent 0.3%</li> </ul> <p>Note: (*) Petroleum Free – Manufactured from environmental friendly Cellulose material made from wood pulp.</p>

<b><u>PRODUCT CHARACTERISTICS</u></b>	<p>Excellent bacterial barrier to both airborne and water-borne contamination</p> <ul style="list-style-type: none"> <li>• Outstanding liquid repellency – water, alcohol, iodine</li> <li>• Exceptional tensile strength – wet or dry</li> <li>• Even texture to allow immediate visual detection of small pinholes or tears</li> <li>• Meets FAA (Flammable Fabrics Act) requirements</li> </ul> <p>Comply with EN ISO 11607-Part 1 and EN 868 – Part 2</p>
---------------------------------------	---

<b><u>PRODUCT USES</u></b>	<ul style="list-style-type: none"> <li>• All types of sterilization center and surgical unit packaging (validated for all steam sterilization methods, including both prevac and gravity equipment)</li> <li>• Inner and outer wrap for medium and large size trays and midweight packs</li> <li>• CSR Wraps are single use item</li> <li>• Use with EtO, steam gamma irradiation, E-beam sterilization technology</li> </ul>
----------------------------	---

<b><u>PACKAGING REQUIREMENTS</u></b>	<table border="0" style="width: 100%;"> <tr> <td style="width: 50%;">Item # SW12 – 12" x 12" – 1,000/Case</td> <td style="width: 50%;">Item # SW15 – 15" x 15" – 500/Case</td> </tr> <tr> <td>Item # SW18 – 18" x 18" – 500/Case</td> <td>Item # SW20 – 20" x 20" – 500/Case</td> </tr> <tr> <td>Item # SW24 – 24" x 24" – 500/Case</td> <td>Item # SW30 – 30" x 30" – 250/Case</td> </tr> <tr> <td>Item # SW36 – 36" x 36" – 250/Case</td> <td>Item # SW40 – 40" x 40" – 250/Case</td> </tr> <tr> <td>Item # SW 45 – 45" x 45" – 250/Case</td> <td>Item # SW48 – 48" x 48" – 100/Case</td> </tr> <tr> <td>Item # SW54 – 54" x 54" – 100/Case</td> <td></td> </tr> </table>	Item # SW12 – 12" x 12" – 1,000/Case	Item # SW15 – 15" x 15" – 500/Case	Item # SW18 – 18" x 18" – 500/Case	Item # SW20 – 20" x 20" – 500/Case	Item # SW24 – 24" x 24" – 500/Case	Item # SW30 – 30" x 30" – 250/Case	Item # SW36 – 36" x 36" – 250/Case	Item # SW40 – 40" x 40" – 250/Case	Item # SW 45 – 45" x 45" – 250/Case	Item # SW48 – 48" x 48" – 100/Case	Item # SW54 – 54" x 54" – 100/Case	
Item # SW12 – 12" x 12" – 1,000/Case	Item # SW15 – 15" x 15" – 500/Case												
Item # SW18 – 18" x 18" – 500/Case	Item # SW20 – 20" x 20" – 500/Case												
Item # SW24 – 24" x 24" – 500/Case	Item # SW30 – 30" x 30" – 250/Case												
Item # SW36 – 36" x 36" – 250/Case	Item # SW40 – 40" x 40" – 250/Case												
Item # SW 45 – 45" x 45" – 250/Case	Item # SW48 – 48" x 48" – 100/Case												
Item # SW54 – 54" x 54" – 100/Case													

<b><u>STORAGE REQUIREMENTS</u></b>	Store in a cool, dry location always from strong sunlight or UV light (30 to 60% RH; 50°F to 104°F / 10°C to 40°C)
------------------------------------	--

<b><u>SHELF LIFE</u></b>	When properly sealed, sterility is maintained up to (5) five years provided the integrity of the roll has not been compromised.
--------------------------	---