



SINGLE USE DISPOSABLES



Articulating and Occlusion Testing Products

- Combined use of articulating film, foil and paper to mark high and low occlusion points provides the most accurate way to mark every point of articulation
- Proprietary manufacturing process combined with stringent quality control guarantees a clear, positive marking without smudging of wet or dry fields
- Color-coded packaging and product descriptions make identification easy
- Articulating paper available in straight, curved, curved with tabs and horseshoe
- Articulating film and foil available in straight

Ordering Information

ARTICULATING PAPER

12 sheets per book/12 books per box

Description	Ref. #
Extra Thin, Blue 0.0015"/38 µm	TPXT
XX Thin, Blue 0.00125"/32 µm	TPXXT
Thick Blue 0.005"/127 µm	TPTH
Thin Blue 0.0028"/71 µm	TPT
Blue Red Combo, 0.0028"/71 µm	TPBR
Blue Red Combo, Thin, Curved 0.003"/71 µm	TPCBRT
Blue Red Combo, with tab, Curved 0.004"/89 µm	TPCBR
Horseshoe Articulating Paper, 12 sheets per book/6 books per box	TPH
Blue Red Combo 0.0035"/89 µm	

ARTICULATING FILM

25 sheets per book/12 books per box

Description	Ref. #
Black Red Combo, Precisely True 0.0005"/13µm	TFMBKR
Blue Thin Precisely True, 0.0005"/13µm	TPP
Red Precisely True, 0.0005"/13µm	TFMR
Blue Red Combo, Precisely True 0.0005"/13µm	TFMBR

ARTICULATING FOIL


25 sheets per book/12 books per box

Description	Ref. #
Blue Red Combo, 0.0005"/12 µm	TFOBR
Red, 2 side 0.0005"/12 µm	TFOR2
Blue, 2 side 0.0005"/12 µm	TFOB2
Blue, 1 side 0.0005"/12 µm	TFOB

NOT MADE WITH NATURAL RUBBER LATEX



Articulating paper available in straight, curved, curved with tab and horseshoe shaped.



Articulating film available in single and double sided straight formats.



Articulating foil available in single and double sided straight formats.

