



# SINGLE-USE DISPOSABLES



## UltraPure® Non-Woven Sponges 4 Ply, Non-Sterile

*A great way to reduce your carbon footprint!*

Eco-friendly UltraPure Non-Woven Sponges. They're soft, lint-free and absorb faster than synthetics, rayons and cotton blends. And, unlike typical synthetic non-woven sponges, they're biodegradable. So you can provide improved treatment for your patients...and protect the environment!

- Made from 100% USA-grown cotton, providing increased comfort and faster absorbency when compared to rayon/polyester and cotton blends
- Hypoallergenic, no toxic chemicals used in manufacturing process; all natural whitening process (H<sub>2</sub>O<sub>2</sub> whitening) without use of bleach and/or chlorine
- Exclusive automated manufacturing process reduces cross-contamination risks, providing the lowest bioburden levels
- Ideal for placement on incision/extraction sites as it easily releases, for increased patient comfort
- Lint-free – extremely soft, smooth, strong and will not adhere to instruments
- Biodegradable – environmentally friendly
- Available in 2 sizes

### Ordering Information

Description	Quantity	Ref. #
2" x 2" / 5 x 5 cm, 45 GR WT*	5000/case	ENC100
4" x 4" / 10 x 10 cm, 40 GR WT*	2000/case	ENC4100



Not made with natural rubber latex

## Ultra Gauze® Sponges 4 Ply, Non-Sterile

- More absorbent than traditional sponges, exclusive 70/30 polyester formed non-woven material replicates the look of woven gauze
- Ideal for placement on incision/extraction sites as it easily releases, for increased patient comfort
- Extremely soft, smooth, strong—virtually lint-free
- Ideal for sustained patient use, disinfection of surfaces, general cleaning
- Available in multiple sizes of 40 gram weight\*

### Ordering Information

Description	Quantity	Ref. #
2" x 2" / 5 x 5 cm	5000/case	ENCNWU
2" x 2" / 5 x 5 cm	8000/case	ENCNWUB
3" x 3" / 7.62 x 7.62 cm	4000/case	ENC3NWU
4" x 4" / 10 x 10 cm	2000/case	ENC4NWU



Not made with natural rubber latex

\* Gram Weight: Weight per square meter of rayon/poly fabric.





## Non-Woven Sponges 4 Ply, Non-Sterile

- Unique rayon/polyester blend
- Provides maximum absorbency
- Less adhesion to wounds for faster healing
- Soft, smooth, strong, virtually lint-free
- Ideal for sustained patient use, disinfecting, general cleaning
- Available in Premium, Advantage Plus and Advantage in multiple gram weights and sizes

### Ordering Information

Description	Quantity	Ref. #
<b>PREMIUM™ 4 Ply Non-Sterile</b> <small>MADE IN USA</small>		
2" x 2" / 5 x 5 cm – 50 GR WT*	5,000/Case	ENCNW
3" x 3" / 7.5 x 7.5 cm – 40 GR WT*	4,000/Case	ENC3NW
2" x 2" / 5 x 5 cm – 50 GR WT*	4,000/Case	ENCNW4K
4" x 4" / 10 x 10 cm – 40 GR WT*	2,000/Case	ENC4NW
<b>ADVANTAGE PLUS® 4 Ply Non-Sterile</b> <small>MADE IN USA</small>		
2" x 2" / 5 x 5 cm – 35 GR WT*	5,000/Case	ENCNWL
<b>ADVANTAGE™ 4 Ply Non-Sterile</b>		
2" x 2" / 5 x 5 cm – 30 GR WT*	5,000/Case	ENCNWLA
3" x 3" / 7.5 x 7.5 cm – 30 GR WT*	4,000/Case	ENC3NWLA
4" x 4" / 10 x 10 cm – 30 GR WT*	2,000/Case	ENC4NWLA

\* Gram Weight: Weight per square meter of rayon/poly fabric.



Not made with natural rubber latex

## Cotton-Filled Exodontia Sponges - 4 Ply

- Made from USP cotton gauze
- Pure, virgin white cotton with 100% cotton-fiber filling
- Highly absorbent
- Available in Premium Non-Sterile, Premium Sterile and Advantage Plus Non-Sterile



Not made with natural rubber latex

### Ordering Information

Description	Quantity	Ref. #
<b>PREMIUM™ Non-Sterile</b>		
2" x 2" / 5 x 5 cm - 125 GR WT*	5,000/Case	ENCC
3" x 3" / 7.5 x 7.5 cm - 125 GR WT*	4,000/Case	ENCC3
4" x 4" / 10 x 10 cm - 125 GR WT*	2,000/Case	ENCC4
<b>PREMIUM™ Sterile</b>		
2" x 2" / 5 x 5 cm - 125 GR WT*	5000/case	ESCC
<b>ADVANTAGE PLUS® Non-Sterile</b>		
2" x 2" / 5 x 5 cm - 90 GR WT*	5,000/Case	ENACC

\* Gram Weight: Weight per sleeve.

## All-Gauze Sponges Premium, Non-Sterile

- Made from USP 100% gauze
- Available in 8 ply, 12 ply or 16 ply
- Highly absorbent

### Ordering Information

Description	Quantity	Ref. #
2" x 2" / 5 x 5 cm - 12 ply	8,000/Case	ENC212
2" x 2" / 5 x 5 cm - 8 ply	5,000/Case	ENC2
3" x 3" / 7.5 x 7.5 cm - 12 ply	4,000/Case	ENC312
3" x 3" / 7.5 x 7.5 cm - 8 ply	4,000/Case	ENC38
4" x 4" / 10 x 10 cm - 16 ply	2,000/Case	ENC416
4" x 4" / 10 x 10 cm - 12 ply	2,000/Case	ENC412
4" x 4" / 10 x 10 cm - 8 ply	4,000/Case	ENC48



Not made with natural rubber latex